


# Применение МУРН для исследования морфологии и внутренней структуры органогелей



Yulia Gorshkova


*FLNP JINR, Dubna, Russia*



MAGDALENA ORDON,  
MIROSŁAWA OSSOWSKA-CHRUŚCIEL  
*Siedlce University of Natural Sciences and  
Humanities, Institute of Chemistry, Siedlce, Poland*

# SMALL ANGLE NEUTRON SCATTERING

## YuMO spectrometer




**q-range**

$7 \times 10^{-3} - 0.5 \text{ \AA}^{-1}$

**size range of object**


**500 – 10 Å (for estimation)**  
**200 Å (max radii of gyration)**



$$\text{wave vector} \quad K_0 = 2 \pi / \lambda$$

$$\text{scattering vector} \quad q = K_0 - k_1$$

- 1 – IBR-2 reactor;
- 2 – first collimator;
- 3 – second collimator;
- 4 – sample;
- 5 – Vn-standard;
- 6 – ring-wire detector;
- 7 – direct beam detector.



**small-angle scattering curve  
(general view)**

## CHALENGE OF INVESTIGATION

### Creation new **SUPRAMOLECULAR ORGANOGENS**

gels with a hydrocarbon dispersion medium

*smart functional*

***nanoscale materials*** with

high potential for a wide range  
of ***advanced applications***

drug delivery


tissue  
engineering

light-  
harvesting  
systems


optoelectronics

# 7OPhOLCA: SCHEME OF SYNTHESIS

**lithocholic acid**




$\xrightarrow{\text{CH}_3\text{OH}}$




**methyl ester of lithocholic acid**

1.  $\text{C}_7\text{H}_{15}\text{O}-\text{C}_6\text{H}_4-\text{OH}$   
2. water/methanol (1:1)



**4-heptyloxyphenyl-lithocholic acid ether (Acronym 7OPhOLCA)**





easy form a hydrogen bonds with the proton acceptor molecules

# GEL FORMATION

	<b>1</b>	<b>2</b>
Chloroform	S	S
Mesitylene	S	P
1,2-Dichlorobenzene	TG (CGC 2 mM)	TG (W)
Chlorobenzene	TG (CGC 2 mM)	GP
Benzene	I	P
Toluene	I	GP
Isopropanol	S	S
<b>DMSO/water</b>	P	<b>TG</b> (CGC 5 mM)
<b>DMF/water</b>	S	<b>TG</b>
MeOH/water	P	S
Acetone/water	S	GP
AcOH/water	S	GP
Dioxane/water	S	TLG
CH <sub>3</sub> CN/water	S	OG

for two compounds:  
**1 – ionic; 2 – non ionic**



**TG** - transparent gel; TLG - translucent gel; GP - gelatinous precipitate;  
S - solution; I - insoluble; P - precipitate; OG - opaque gel; W - weak

do not form gel with bile acids the solvents with high molecular weight or  
without proton donor-acceptor groups

form gel in mixtures of DMSO (dimethyl sulfoxide)/water (1:2 to 3:2),  
DMF (dimethylformamide)/water (2:3 to 3:2), 1,4-dioxane/water (1:4)

[Maitra, J. Org. Chem. 2011].

*DMSO – organic solvent have proton acceptor group.  
Gel systems are formed by hydrogen bonds.*

# SULFOXIDES: PROPERTIES AND APPLICATIONS

## CRYOPROTECTOR

### PROPERTY

The ability to bind free radicals



Reducing damage of biological objects

### APPLICATION

Cell biology  
Cryobiology  
Medicine  
Agriculture

## DMSO


## SOLVENT

### PROPERTY

High permeability through biological membranes

### APPLICATION

Medicine  
Pharmacology



hydrophobic interactions play a crucial role in the interaction

# SCIENTIFIC INTEREST


## Morphology

## Internal structure


## Phase transitions



Solutions of gel 7OPhOLCA in DMSO- $d_6$  in concentrations: 1 – 0.005g/ml; 2 – 0.010g/ml; 3 – 0.015g/ml; 4 – 0.020g/ml; 5 – 0.025g/ml at room temperature



# DSC METHOD




Thermograms of full cycle heating and cooling for 70PhOLCA in DMSO-d<sub>6</sub> (0.015g/ml), in the range of -20°C to 60°C with rate of 5deg/min for sample with mass 14.02mg

# SANS METHOD: Effect of concentration

Guinier region

Morphology

Internal structure




The scattering intensity is determined by a power law:  $I(q) \sim q^{-m}$ , where  $m$  is a power-law exponent related to a fractal dimension as  $D = 6 - m$ . The gradient gives information about the morphology of the sample.

The scattering on the individual lamellas with  $d = 43.4 \pm 0.1 \text{\AA}$

SANS curves for 70PhOLCA in DMSO- $d_6$  at concentrations  $C = 0.005, 0.010, 0.015, 0.020$  and  $0.025 \text{ g/ml}$  at  $T = 10^\circ\text{C}$

# SANS METHOD: Temperature effect



$m = -4.0 \pm 0.2$  – objects with a smooth surface,


$m = -3.6 \pm 0.1$  – surface fractal,

$m = -2.0 \pm 0.1$  – a system consist from the fully separated lamellas



SANS curves for 70PhOLCA in DMSO-d<sub>6</sub> with C = 0.015g/ml at T = 10, 25 and 45°C

*Increasing of the temperature led to the morphology transitions  
Gel 1 → Gel 2 and Gel 2 → Gel 3.*

# DSC METHOD: Morphology transition in Gel phase




# SANS METHOD: Sol – Gel transition



Proposed Sol – Gel transition scheme for the 7OPhOLCA in DMSO-d<sub>6</sub> determined by SANS. The values 4, 3.0 - 3.9, 2 and 1 correspond to the objects with a smooth surface, surface fractal, lamellas in solution and disordered state, respectively

Gel – Sol transition is reversible with hysteresis on the temperature



# PHASE TRANSITION OF PC MEMBRANES IN THE SULFOXIDES PRESENTS



Phase diagrams of the **spontaneous**  
MLVs DMPC (2 wt %)


in DMSO/D<sub>2</sub>O  
 $X_{\text{DMSO}} = 0.2$

in DESO/D<sub>2</sub>O  
 $X_{\text{DESO}} = 0.2$




Transition  
 $L_{\beta'} \rightarrow L_\alpha$


Transitions  
 $L_{\beta'} \rightarrow P_{\beta'} \rightarrow L_\alpha$



# SANS METHOD: Time effect



**A**  
Optical microscopy image for 70PhOLCA in DMSO-d<sub>6</sub> at concentration  
C = 0.025 g/ml (**A**) after preparation, (**B**) after 4 days with max  
magnification 6.3x of the crystallization center and (**C**) after 14 days



SANS curves for 70PhOLCA in DMSO-d<sub>6</sub> with C = 0.015, 0.020 and  
0.025g/ml at T=20°C after preparation (**left**) and after 14 days (**right**)

# INFRARED SPECTROSCOPY (FT - IR)




Figure of full-width at half-maximum (FWHM) in function of temperature of hydroxyl group, during heating. A – stretching vibrations ( $\nu_{OH}$ ); B – deformation vibrations ( $\delta_{OH}$ ) in the range of temperature from 30  $^{\circ}\text{C}$  to 60  $^{\circ}\text{C}$  for sample "3".

# SUMMARY

Type of transition	HEATING CYCLE						
	0.015g/ml			0.025g/ml			
	DSC	SANS	FTIR	DSC	SANS		
	5deg/min			2deg/min			
	14.02mg			26.21mg			
Cr - Gel	T <sub>p</sub>	4.72	X	X	4.37	X	
	ΔH	61.21			52.99		
	ΔS	0.22028			0.19094		
Gel - Fractal form of gel	T <sub>p</sub>	32.92	X	X	X	25	
	ΔH	0.20			X		
	ΔS	0.00065			X		
Fractal form of gel - Lamellar form of gel	T <sub>p</sub>	42.25	45	42	X	28	
	ΔH	0.03			X		
	ΔS	0.00010			X		
Lamellar form of gel - Sol	T <sub>p</sub>	X	52	50	33.60	35	
	ΔH	X			1.29		
	ΔS	X			0.00421		
COOLING CYCLE							
Type of transition	0.015g/ml			0.025g/ml			
	DSC	SANS	FTIR	DSC	SANS		
	5deg/min			2deg/min			
	14.02mg			26.21mg			
	T <sub>p</sub>	18.58	18	X	X	29	
Sol - Lamellar form of gel	ΔH	0.21			X		
	ΔS	0.00072			X		
Lamellar form of gel – Fractal form of gel	T <sub>p</sub>	14.99	15	X	X	21	
	ΔH	0.08			X		
	ΔS	0.00028			X		
Fractal form of gel - Gel	T <sub>p</sub>	11.46	11	X	15.29	18	
	ΔH	0.09			1.18		
	ΔS	0.000316			0.00409		
Gel – Cr	T <sub>p</sub>	-6.41	X	X	-4.95	X	
	ΔH	81.90			69.64		
	ΔS	0.29296			0.25041		

# USER POLICY IN FLNP JINR

## IBR-2 Spectrometers:

**Diffraction:** HRFD, DN-2, DN-12, DN-6, SKAT-EPSILON, FSD



**Small-angle scattering:** YuMO

**Reflectometry:** REMUR, REFLEX, GRAINS

**Inelastic scattering:** DIN-2PI, NERA

**Nuclear Physics:** ISOMER, KOLHIDA

**Irradiation Facility**



Period for proposal submission  
**September 1 - October 15**  
March 1 - April 15

<http://ibr-2.jinr.ru>

# Conferences on Condensed Matter Research at IBR-2 reactor

**International Conference**

**Condensed Matter Research at the IBR-2**

**October 11 - 15, 2015**

**Dubna, Russia**

**TOPICS**

- Functional and nanostructured materials;
- Magnetic colloid systems;
- Layered magnetic nanostructures;
- Carbon nanostructures;
- Materials under extreme conditions;
- Soft condensed matter (biological nanosystems, lipid membranes, polymers);
- Lattice and molecular dynamics of materials;
- Texture and properties of rocks and minerals;
- Residual stresses in materials and products;
- Neutron imaging;
- Development of IBR-2 instruments;
- Development of neutron scattering techniques;
- Development of neutron detectors

**ORGANIZING COMMITTEE**

Prof. D.P. Kostenko – Chairman  
Dr. M.V. Avdeev  
Prof. A.M. Balagurov  
Dr. G.D. Bokuchava  
Dr. D.M. Chudoba  
Mrs. T.S. Donskova – Secretary  
Mrs. Yu.E. Gorishkova – Scientific Secretary  
Mr. S.V. Kosenkov  
Dr. B.N. Savchenko

**PROGRAMME COMMITTEE**

Prof. V.L. Alsanov (Russia)  
Prof. P.A. Aleksiev (Russia)  
Prof. A.M. Balagurov (Russia)  
Prof. P. Balogý (Slovak Republic)  
Prof. A.V. Belushkin (Russia)  
Prof. E. Burzo (Romania)  
Prof. H. Fuss (Germany)  
Prof. V.G. Kanter (Moldova)  
Prof. D. Nagy (Hungary)  
Prof. W. Nawrocki (Poland)  
Prof. V.A. Somerkov (Russia)

**DEADLINES**

- First announcement: 15.04.2015
- On-line registration and abstract submission: 30.06.2015
- Visa support: 30.06.2015
- Programme: 20.09.2015
- Arrival: 11.10.2015
- Departure: 15.10.2015

**CMR @ IBR-2**

E-mail: [cmr@nf.jinr.ru](mailto:cmr@nf.jinr.ru)  
Web Site: <http://cmr-ibr.jinr.ru>  
Tel.: (+7-496-21)62774  
Fax: (+7-496-21) 65484

*Thank  
you for  
your  
attention!*



*Up to  
new  
meetings  
in Dubna*



E-mail: [cmr@nf.jinr.ru](mailto:cmr@nf.jinr.ru)  
Web Site: <http://cmr-ibr.jinr.ru>  
Fax: (+7-496-21) 65484



Yu. Gorshkova  
FLNP JINR

# WHAT WE CAN LEARN USING SANS: OBJECTS AND PARAMETERS

*shape and size  
of particles*

*structure*

*volume  
fraction*

*size  
distribution*

*molecular weight*


*thickness of the layers  
and repeat spacing*

*particles  
density*

*mosaicity*

*fractal  
dimension*

*aggregation number*



## Main parameters of YuMO instrument

Parameters	Value
Flux on the sample (thermal neutrons)	$10^7 - 4 \times 10^7 \text{ n}/(\text{s cm}^2)$ [1]
Used wavelength	0.5 Å to 8 Å #
Q-range	$7 \times 10^{-3} - 0.5 \text{ \AA}^{-1}$
Dynamic Q-range	$q_{\max}/q_{\min}$ up to 100
Specific features	Two detectors system [2,3], central hole detectors
Size range of object *	500 – 10 Å
Intensity (absolute units -minimal levels)	0.01 cm <sup>-1</sup>
Calibration standard	Vanadium during the experiment [4]
Size of beam on the sample	8 – 22 mm <sup>2</sup> @
Collimation system	Axial
Detectors	He <sup>3</sup> -fulfilled, home made preparation, 8 independent wires [5]
Detector (direct beam)	<sup>6</sup> Li-convertor (home made preparation)
Condition of sample	In special box in air
Q-resolution	low, 5-20%
Temperature range	-50°C -+130°C ^ (Lauda)
Temperature range	700 °C ** (Evrotherm)
Number of computer controlled samples	14 ***
Background level	0.03 – 0.2 cm <sup>-1</sup>
Mean time of measurements for one sample	1 h +
Frequency of pulse repetition	5 Hz
Electronic system	VME
The instrument control software complex	SONIX [6]

# WHAT WE CAN LEARN USING SANS: SIZE RANGE

